



中国认可
国际互认
检测
TESTING
CNAS L0095

Page 1 of 8 Pages

No.: RZUN2023-6742-M1

检测报告

TEST REPORT

UN38.3

NAME OF SAMPLE:

Li-MnO₂ Button Cell

产品名称:

锂一二氧化锰扣式电池

CLIENT:

CHANGZHOU JINTAN CHAOCHUANG BATTERY
CO.,LTD

委托单位:

常州金坛超创电池有限公司

CLASSIFICATION OF TEST:

Commission Test

检测类别:

委托测试

威凯检测技术有限公司
CVC Testing Technology Co., Ltd.

检测报告

TEST REPORT

No.: RZUN2023-6742-M1

Page 2 of 8 Pages

Name of samples: Li-MnO ₂ Button Cell 样品名称: 锂-二氧化锰扣式电池	Type/Model: 型号规格: CR2016 3V 75mAh
Color: Silver 样品颜色: 银色	Physical shape: Button 样品形状: 纽扣
Commissioned by: CHANGZHOU JINTAN CHAOCHUANG BATTERY CO.,LTD 委托单位: 常州金坛超创电池有限公司	Commissioner address: XIYANG VILLAGE MAOLU TOWN JINTAN CHANGZHOU JIANGSU 委托单位地址: 江苏省金坛区茅麓镇西阳村
Manufacturer: CHANGZHOU JINTAN CHAOCHUANG BATTERY CO.,LTD 制造商: 常州金坛超创电池有限公司	Manufacturer address: XIYANG VILLAGE MAOLU TOWN JINTAN CHANGZHOU JIANGSU 制造商地址: 江苏省金坛区茅麓镇西阳村
Factory: CHANGZHOU JINTAN CHAOCHUANG BATTERY CO.,LTD 生产厂: 常州金坛超创电池有限公司	Factory address: XIYANG VILLAGE MAOLU TOWN JINTAN CHANGZHOU JIANGSU 生产厂地址: 江苏省金坛区茅麓镇西阳村
Classification of test: Commission Test 检测类别: 委托测试	Quantity of sample: 1 cell 样品数量: 1 个电芯
Tested according to: 测试标准: ST/SG/AC.10/11/Rev.7/Amend.1/Section 38.3	Sample identification: 样品标识序号: 1#
Receiving date: 接样日期: 2024-12-26	Means of receiving: Submitted by commissioner 接样方式: 委托单位送样
Completing date: 完成日期: 2025-01-02	Test item: - 测试项目: -
Test conclusion: 检测结论: The Li-MnO₂ Button Cell is submitted by CHANGZHOU JINTAN CHAOCHUANG BATTERY CO.,LTD. The sample's differences between this report and original report RZUN2023-6742 are comply with the requirements of the Section 38.3 of the Seventh revised edition Amendment 1 of the Manual of Tests and Criteria (ST/SG/AC.10/11/Rev.7/Amend.1/Section 38.3). No need to test. 由常州金坛超创电池有限公司送检的锂-二氧化锰扣式电池, 与原报告 RZUN2023-6742 的样品差异符合联合国《试验和标准手册》第七修订版修正 1 第 38.3 节相关规定, 无需测试。 Seal of CVC CVC 盖章 Date of issue: 签发日期: 批准签发日期	

Title: Manager
批准人职务: 经理

Approved by: Zhang Siyao Reviewed by: Zheng Feida Tested by: Jiang Weiwei

批 准: 批准签名 EN 审 核: 审核签名 EN 检 测: 编制签名 EN

Description and illustration of the sample:

样品说明及描述:

The sample's status is good

样品状况良好。

Aggregate Lithium content of each cell/ 单个电芯的锂含量: 0,025g

Description of the sampling procedure:

取样程序的说明:

/

Description of the deviation from the standard, if any:

试验结果不符合标准项的说明:

/

Remarks:

备注:

1. Throughout this report a comma is used as the decimal separator.
本报告中以逗号代替小数点。
2. This report is transformed from original report RZUN2023-6742, which is issued on 2023.12.28, no new test item is involved.
本报告是 RZUN2023-6742 的变更报告, 原报告签发于 2023 年 12 月 28 日, 不涉及新的测试项目。
3. New trademarks have been added to this report, the modified items did not change the key materials, design and production process of the sample of original report, the production factory entity has no change.
本报告中增加了新商标, 该变更未改变原报告中样品的关键材料, 产品设计和生产工艺, 生产厂也未改变。
4. This report needs to be used with the original report.
本报告需和原报告一起使用。
5. Trade mark/商标: CR2016 (无品牌)、CR2016 (GNENG)、CR2016 (STAR BULLY)、CR2016 { (CHAOCHUANG) (SC100、SC200、SC360、SC500、SC800、SC628、SC508、SC928、SC358、SC218、SC158)、CR2016UL 桶、CR2016mAh 桶、CR2016 老中、CR2016STD、CR2016 老中桶、CR2016 中性桶、CR2016 大桶、CR2016 CHAOCHUANG、CR2016 S CHAOCHUANG、CR2016 L CHAOCHUANG、CR2016 H CHAOCHUANG、CR2016 H ADMA、CR2016 CHAOCHUANG LED、CR2016 S ADMA、CR2016 ADMA、CR2016 H、CR2016 TOSHIBA、CR2016 REANTA、CR2016BEXEL、CR2016 BIXLINK、CR2016 老中 2、CR2016 MAXELL、CR2016 MORITOKU、CR2016 Navigator、CR2016 双鹿中性、CR2016 金力中性、CR2016 SUNNYBATT、CR2016 GoldenPower、CR2016 南孚、CR2016 GIGAMAX、CR2016 ANSMANN、CR2016 MastAK、CR2016 HW、CR2016 TDM、CR2016 KEYKO、CR2016 艾科夫、CR2016 ERGOLUX、CR2016 Focusray、CR2016 HT、CR2016 RAYOVAC、CR2016 Westinghouse、CR2016 CamelionCR2016 Excell、CR2016 perfeo、CR2016 UNOMRT、CR2016 FOXLUX、CR2016 TCL、CR2016 CHANGHONG、CR2016 TAURUS、CR2016 HYUNDAI、CR2016 KRISBOW、CR2016 CECASA、CR2016 SONNEN、CR2016 NOMA CR2016 JUPITER、CR2016 PHILIPS、DL/CR2016 DURACELL、CR2016 Starmax、CR2016 Jauch、NOMA、Vomiko、双鹿中、MWK、ZEUS、Ecobattery、KPENTON、intelbras、amazon、PC、金力中、GoldenPower、M&G、Energizer、TGL、MORITOKU、TAURUS、renata、PC、RAYOVAC、amazonbasics、Mirex、TESLA、KOOMAX、PROCELL、DURACELL、HIM、CELEWELL、Camelion 儿安、Matsushima、中性儿安、DALICELL、DOMSEM、M 桶儿安、H 桶儿安、L 桶儿安、老中儿安、三菱、PHILIPS、JAUCH 儿安、Intenso、KODAK、SENCOR、全英文儿安、WILKINSON、KNOWPIANO、H MAXFORT、Heybright、UP&UP、DUNBA、H 超创 CE 桶儿安、超创 CE 桶儿安、CE 桶儿安(CELL)、GutAlkali、DELIPOW、WHITECRANE 儿安、REXANT、KODAK-CE、帕德、Greenlight、ZCELL、LICB。

Photos of Samples and Labels/样品照片及标识

Cell/电芯 (CR2016 3V 75mAh)



Photos of Samples and Labels/样品照片及标识

Cell/电芯 (CR2016 3V 75mAh)



Photos of Samples and Labels/样品照片及标识

Cell/电芯 (CR2016 3V 75mAh)



注 意 事 项 Important

1. 报告无检测单位印章无效。
The test report is invalid without the seal of CVC.
2. 未经本试验室书面同意，不得部分地复制本报告。
Nobody is allowed to photocopy or partly photocopy this test report without written permission of CVC.
3. 本报告无批准人、审核人及检测人签名无效。
The test report is invalid without the signatures of Ratifier, Reviewer and Testing engineer.
4. 本报告涂改无效。
The test report is invalid if altered,
5. 对检测报告若有异议，应于收到报告之日起十五天内向检测单位提出。
Objections to the test report must be submitted to CVC within 15 days.
6. 本报告仅对送检样品负责。
The test report is valid for the tested samples only.
7. 判定栏中“-”表示“不需要判定”，“P”表示“通过”，“F”表示“不通过”，“N/A”表示“不适用”。
As for the Verdict, “-” means “no need for judgement”, “P” means “pass”, “F” means “fail” and “N/A” means “not applicable”.

****报告中未加 CMA 标志时，检测数据和结果仅供科研、教学或内部质量控制之用。****
The test data and test results given in this test report should only be used for purposes of scientific research, teaching and internal quality control when the CMA symbol is not presented.

地 址： 广东省广州市科学城开泰大道天泰一路 3 号

广东省广州市黄埔区光谱东路 179 号百事高智慧园 D 栋（测试地点）

Lab Address: No.3, Tiantai 1st Road, Kaitai Avenue, Science City, Guangzhou, Guangdong, China.

Testing Location: Building D, BASIGO INTELLIGENT, No.179, Guangpu East Road, Huangpu District, Guangzhou, P. R. China.

电 话(Tel): 020 32293888

邮政编码(Post Code): 510663

传 真(FAX): 020 32293889

E-mail: office@cvc.org.cn

<http://www.cvc.org.cn>